

MR Pelvis (Routine Body)

Last Updated: 12/2025

Name	Plane	Sequence	FOV
COR T2 Single Shot	Coronal	T2 Single Shot (AKA: SS TSE / HASTE)	
SAG T2	Sagittal	T2 FSE (AKA: TSE)	28-30 cm (pelvic sidewall to sidewall)
COR T2	Coronal	T2 FSE (AKA: TSE)	28-30 cm
AX T2	Axial	T2 FSE (AKA: TSE)	28-30 cm
AX DWI	Axial	DWI	28-30 cm
AX T1	Axial	DUAL ECHO IN/OUT PHASE GRE	28-30 cm
AX T2 FATSAT	Axial	T2 FAT SAT FSE (AKA: TSE / HASTE)	28-30 cm
SAG DYNAMIC GAD	Axial	LAVA DYNAMIC (INCLUDE PRE-CONTRAST) (AKA: THRIVE / VIBE)	28-30 cm
AX GAD	Sagittal and Coronal	LAVA DELY (ASK: THRIVE / VIBE)	28-30 cm

Notes:

- LAVA sequence is set on multiphase 6 sets. Run the first one before contrast injection.
- Pre, art, late art, portal venous phase, equilibrium 3min, and 5 min delay.
- Test bolus: position at aortic bifurcation, time from first flash of contrast.