

IV Contrast studies: 2ml/kg 100ml max dose

Routine: Preferred - Power Injector Set To Deliver Bolus In 30-40sec.

(Manual Injection Over 30 Seconds)

Delay Scan Start 70 Sec from Start of Injection or if chest only 45 sec from start of injection

Delay kidney or bladder only at radiologist discretion

RECONS:

Name	Type	Window	Slice	Interval	Scan Range
TOPOGRAM					Top of liver to either iliac crest for abd or pubic symphysis for abdpel
AXIAL ST	MPR	Soft Tissue	3 mm	2 mm	
CORONAL ST	MPR	Soft Tissue	3 mm	2 mm	Skin to skin
SAGITTAL ST	MPR	Soft Tissue	3 mm	2 mm	Skin to skin

CT chest for pectus deformity (Haller index): low dose CT chest, no IV contrast **
during end expiration **