

Name	Thick	Space	Comments
<b>AX T2 STIR</b>	3 mm	1 mm	
<b>AX T2 BLADE</b>	3 mm	1 mm	
<b>AX STIR WS Silicone</b>	3 mm	1 mm	
<b>SAG T2 BLADE</b>	3 mm	1 mm	Separate in to left and right breast series
<b>SAG T2 STIR Silicone</b>	3 mm	1 mm	Separate in to left and right breast series
<b>INJECT CONTRAST</b>			
<b>AX POST T1 3D FS</b>	2 mm	0 mm	Check images before proceeding with dynamics
<b>AX POST T1 3D DYN</b>	2 mm	0 mm	60s, 120s, 180s
<b>AX POST SUB</b>	2 mm	0 mm	60s, 120s, 180s
<b>AX POST SUB MIPS</b>			60s, 120s, 180s

Protocol Updated: 12/2025