

CT Abdomen (Multiphasic Liver)

Last Updated: 12/2025

Name	Plane	Type	Window	Thick	Space	Scan Range
SCOUT						
AXIAL PRE	Axial	MPR	Soft Tissue	5 mm	5 mm	Diaphragm through the liver
Inject contrast - 35 second delay						
AXIAL ARTERIAL	Axial	MPR	Soft Tissue	5 mm	5 mm	Diaphragm through the liver
Venous can at 70 seconds after injection						
AXIAL VENOUS	Axial	MPR	Soft Tissue	5 mm	5 mm	Diaphragm to the iliac crests
Delay scan at 5 minutes after injection						
AXIAL DELAY	Axial	MPR	Soft Tissue	5 mm	5 mm	Diaphragm through the liver
COR MPR x 4	Coronal	MPR	Soft Tissue	2.5 mm	2.5 mm	
SAG MPR x 4	Sagittal	MPR	Soft Tissue	2.5 mm	2.5 mm	

Notes:

- Contrast: Per site, 80-100 cc.
- No oral contrast unless specified by the ordering provider.
- Craniocaudal scan direction.