

CT Abdomen Pelvis (Adrenal Mass)

Last Updated: 9/3/2025

Name	Plane	Type	Window	Thick	Space	Scan Range
SCOUT						
AXIAL PRE	Axial	MPR	Soft Tissue	2.5 mm	2.5 mm	Pre: diaphragm to mid kidney
Inject contrast - 70 second delay						
AXIAL POST 70 second delay	Axial	MPR	Soft Tissue	2.5 mm	2.5 mm	Post: diaphragm to iliac crests
15 minute delay						
AXIAL DELAY 15 min	Axial	MPR	Soft Tissue	2.5 mm	2.5 mm	Delay: Diaphragm to mid kidney
COR MPR x 3 phases	Coronal	MPR	Soft Tissue	2.5 mm	2.5 mm	
SAG MPR x 3 phases	Sagittal	MPR	Soft Tissue	2.5 mm	2.5 mm	

Notes:

- Contrast: Per site, 100 cc.
- No oral contrast.
- Craniocaudal scan direction.