

CT Abdomen Pelvis (CT IVP)

Last Updated: 12/2025

Name	Plane	Type	Window	Thick	Space	Scan Range
SCOUT						
AXIAL PRE	Axial	MPR	Soft Tissue	5 mm	5 mm	diaphragm to trochanters
Double bolus technique:						
- 50 cc at time 0						
- 80 cc at 5 min						
- Scan at 7 min from time 0						
AXIAL POST	Axial	MPR	Soft Tissue	5 mm	5 mm	diaphragm to trochanters
COR MPR x 2	Coronal	MPR	Soft Tissue	2.5 mm	2.5 mm	
SAG MPR x 2	Sagittal	MPR	Soft Tissue	2.5 mm	2.5 mm	

Notes:

- Contrast: Per site, 100-120 cc.
- No oral contrast.
- Patients should be instructed to drink approximately 800-1000cc water, as tolerated, 30min prior to the study.