

CT Abdomen Pelvis (Enterography)

Last Update: 12/2025

Name	Plane	Type	Window	Thick	Space	Scan Range
SCOUT						
Inject contrast - 50 second delay						
AX ABD/PEL	Axial	MPR	Soft Tissue	2.5 mm	2.5 mm	Lung base through pubic symphysis
COR MPR	Coronal	MPR	Soft Tissue	2.5 mm	2.5 mm	
SAG MPR	Sagittal	MPR	Soft Tissue	2.5 mm	2.5 mm	

Notes:

- Contrast: Per site, 100 cc.
- Craniocaudal scan direction.
- After a 4 hour fast, the patient is to drink approx. 1800 ml of 0.1% (weight/volume) barium sulfate suspension (Volumen; E-Z Em) divided into four 450ml oral doses given every 20 min, beginning 60 min prior to scanning w/ the last dose when the patient is on the CT table.