

CT Abdomen Pelvis With or Without (Routine)

Last Updated: 12/2025

Name	Plane	Type	Window	Thick	Space	Scan Range
SCOUT						
Inject contrast - 70 second delay						
AX ABD/PEL	Axial	MPR	Soft Tissue	5 mm	5 mm	Lung base through pubic symphysis
COR MPR	Coronal	MPR	Soft Tissue	2.5 mm	2.5 mm	
SAG MPR	Sagittal	MPR	Soft Tissue	2.5 mm	2.5 mm	

Notes:

- Contrast: Per site, 100 cc, 70 second delay.
- Craniocaudal scan direction.
- Oral contrast per institutional protocol or ordering provider's request.