

## CT Chest Super Dimension (Super D)

Last Updated: 9/2025

Name	Plane	Type	Window	Thick	Space	Scan Range
<b>NON-CONTRAST EXAM</b>						
<b>SCOUT</b>						
<b>AX CHEST</b>	Axial	MPR	Soft Tissue	2.5 mm	2.5 mm	2-3 cm above lung apices through adrenal glands
<b>AXIAL DETAIL</b>	Axial	MPR	Lung	1.25 mm	1.25 mm	
<b>AXIAL LUNG</b>	Axial	MPR	Lung	2.5 mm	2.5 mm	
<b>COR MPR</b>	Coronal	MPR	Soft Tissue	2.5 mm	2.5 mm	
<b>SAG MPR</b>	Sagittal	MPR	Soft Tissue	2.5 mm	2.5 mm	

### Notes:

- Contrast: None.
- Indication: Lung mapping prior to bronchoscopy.
- Reference vendor site: <http://www.inlandimaging.com/protocols/ct-chest-super-d>.
- Hyperventilate, then breath-hold in inspiration.
- Craniocaudal scan direction.