

CT Facial Bones

Last Updated: 12/2025

Name	Plane	Type	Window	Thick	Space	Scan Range
SCOUT						AP and Lat Head
AXIAL ST	Axial	MPR	Soft Tissue	2.5 mm	2.5 mm	Mandibular Eminence to Frontal Sinus
AXIAL BONE	Axial	MPR	Bone	2.5 mm	2.5 mm	Mandibular Eminence to Frontal Sinus
CORONAL BONE	Coronal	MPR	Bone	2.5 mm	2.5 mm	Nose to Sella
SAGITTAL BONE	Sagittal	MPR	Bone	2.5 mm	2.5 mm	Skin to Skin

Notes:

- Noncontrast for trauma, contrast for infection or mass
- For CT Face W/ or W/WO- include post contrast 2.5mm axial, coronal and sagittal MPRs in soft tissue window

Contrast:

- Contrast Type: Per site
- Contrast Dose: 75-100 cc
- Delay: 60 sec