

CT Head Stereotactic/Stryker

Last Updated: 12/2025

Name	Plane	Type	Window	Thick	Space	Scan Range
SCOUT						AP and Lat Head
INJECT CONTRAST - 2 MINUTE DELAY						
AXIAL ST	Axial	MPR	Soft Tissue	1 mm	1 mm	Mandible to Vertex
AXIAL BONE	Axial	MPR	Bone	1 mm	1 mm	Mandible to Vertex
CORONAL	Coronal	MPR	Soft Tissue	1 mm	1 mm	Skin to Skin
SAGITTAL	Sagittal	MPR	Soft Tissue	1 mm	1 mm	Skin to Skin

Notes:

- All acquisitions must be helical
- Expanded FOV from mandibular eminence to scalp vertex and including face/nose

Contrast:

- Contrast Type: Per site.
- Contrast Dose: 75-100 mL
- Delay: 2 min