

## CT Head With

Name	Plane	Type	Window	Thick	Space	Scan Range
SCOUT						AP and Lat Head
<b>INJECT CONTRAST - 2 MINUTE DELAY</b>						
AXIAL ST	Axial	MPR	Soft Tissue	5 mm	5 mm	Skull Base to Vertex
AXIAL BONE	Axial	MPR	Bone	5 mm	5 mm	Skull Base to Vertex
CORONAL	Coronal	MPR	Soft Tissue	5 mm	5 mm	Skin to Skin
SAGITTAL	Sagittal	MPR	Soft Tissue	5 mm	5 mm	Skin to Skin

### Notes:

- All acquisitions should be helical unless unavailable

### Contrast:

- Contrast Type: Per site
- Contrast Dose: 75-100 mL
- Delay: 2 min