

CT Parathyroid for UTRGV

Last Updated: 12/2025

| Name | Plane | Type | Window | Thick | Space | Scan Range |
|------------------------------------------------|-------|------|--------|---------|---------|------------|
| Pre-Contrast | | | | | | |
| Axial ST | Ax | MPR | ST | 1.25 mm | 1.25 mm | see below |
| Cor ST | Cor | | | | | |
| Sag ST | Sag | | | | | |
| Inject contrast 25 sec delay (Arterial) | | | | | | |
| Axial ST | Ax | MPR | ST | 1.25 mm | 1.25 mm | see below |
| Cor ST | Cor | | | | | |
| Sag ST | Sag | | | | | |
| 65 sec delay (Venous) | | | | | | |
| Axial ST | Ax | MPR | ST | 1.25 mm | 1.25 mm | see below |
| Cor ST | Cor | | | | | |
| Sag ST | Sag | | | | | |

Notes:

Scan from the bottom of the mandible through the carina

[PARATHYROID-4D-GE-REVOLUTION-EVO-64-SLICE-72-KW](#)