

CT Pelvis (Cystogram)

Last Updated: 12/2025

Name	Plane	Type	Window	Thick	Space	Scan Range
SCOUT						
AXIAL PELVIS PRE	Axial	MPR	Soft Tissue	2.5 mm	2.5 mm	1 cm above the iliac crest through the pubic symphysis
Instill 30 cc of Omnipaque 300 into 500cc of normal saline. Retrograde drip into the bladder, via the patient's foley, either 500 cc or until the patient cannot tolerate.						
Do not ever inject with pressure or squeeze bag as bladder injury could occur or worsen.						
AXIAL PELVIS FILL	Axial	MPR	Soft Tissue	2.5 mm	2.5 mm	
COR MPR x 2	Coronal	MPR	Soft Tissue	2.5 mm	2.5 mm	
SAG MPR x 2	Sagittal	MPR	Soft Tissue	2.5 mm	2.5 mm	

Notes:

- Contrast: Per site.
- Craniocaudal scan direction.
- Instill 30 cc of Omnipaque 300 into 500cc of normal saline. Retrograde drip into the bladder, via the patient's foley, either 500 cc or until the patient cannot tolerate. Do not ever inject with pressure or squeeze bag as bladder injury could occur or worsen.
- No oral contrast.