

**CT Sacrum/Coccyx**

Last Updated: 12/2025

Name	Plane	Type	Window	Thick	Space	Scan Range
SCOUT						AP and Lat
AXIAL ST	Axial	MPR	Soft Tissue	2 mm	2 mm	L5 Through Coccyx
AXIAL BONE	Axial	MPR	Bone	2 mm	2 mm	L5 Through Coccyx
CORONAL BONE	Coronal	MPR	Bone	2 mm	2 mm	Skin to Ant Margin of Sacrum
SAGITTAL BONE	Sagittal	MPR	Soft Tissue	2 mm	2 mm	SI Joint to SI Joint

**Notes:**

- Bone algorithm/kernel
- Noncontrast for trauma and DDD. Contrast for indication of infection or tumor.
- For CT Sacrum W/ or W/WO- include post contrast 3mm axial, coronal and sagittal MPRs in soft tissue windows

**Contrast:**

- Contrast Type: Per site
- Contrast Dose: 75-100 cc
- Delay: 60 sec