

**CT Sinus**

Last Updated: 12/2025

Name	Plane	Type	Window	Thick	Space	Scan Range
SCOUT						AP and Lat Head
AXIAL ST	Axial	MPR	Soft Tissue	2 mm	2 mm	Top of Frontal Sinus to Maxillary Dentition
AXIAL BONE	Axial	MPR	Bone	2 mm	2 mm	Top of Frontal Sinus to Maxillary Dentition
CORONAL BONE	Coronal	MPR	Bone	2 mm	2 mm	Nose to Sella
SAGITTAL BONE	Sagittal	MPR	Bone	2 mm	2 mm	Skin to Skin

**Notes:**

- Orient axial plane parallel to hard palate
- Noncontrast unless indication for mass
- For CT Sinus W/ or W/WO- include post contrast 2mm axial, coronal and sagittal MPRs in soft tissue window

**Contrast:**

- Contrast Type: Per site
- Contrast Dose: 75-100 cc
- Delay: 60 sec