

CT Venogram Head

Last Updated: 12/2025

Name	Plane	Type	Window	Thick	Space	Scan Range
SCOUT						AP and Lat Head
AXIAL ST	Axial	MPR	Soft Tissue	1 mm	1 mm	Scalp vertex to C1
CORONAL ST	Coronal	MPR	Soft Tissue	2 mm	2 mm	Skin to Skin
SAGITTAL ST	Sagittal	MPR	Soft Tissue	2 mm	2 mm	Skin to Skin
CORONAL MIP	Coronal	MIP	Soft Tissue	10 mm	2 mm	Skin to Skin
SAGITTAL MIP	Sagittal	MIP	Soft Tissue	10 mm	2 mm	Skin to Skin

Notes:

Contrast:

- Contrast Type: Per site
- Contrast Dose: 75-100 cc
- Delay: 45 sec