

## CT Whole Body Myeloma Screening

Last Updated: 5/2025

Name	Plane	Window	Thick	Space	Scan Range
AX Bone	AX	BONE	2.5 mm	2.5 mm	
AX Chest	AX	ST	2.5 mm	2.5 mm	Top of skull through knees
COR MPR	AX	BONE	2.5 mm	2.5 mm	
SAG MPR	AX	BONE	2.5 mm	2.5 mm	

### Notes:

- Noncontrast exam
- Use bone algorithm for bone windows
- Low dose protocol: tube voltage: 120 kV, tube current: 40-50 mAs
- Collimation: 16.0 x 1.5