

CTA Abdomen Pelvis (Aortic Stent Graft/Endoleak)

Last Updated: 12/2025

Name	Plane	Type	Window	Thick	Space	Scan Range
SCOUT						
AXIAL PRE	Axial	MPR	Soft Tissue	2 mm	2 mm	Lung base through pubic symphysis
Inject contrast - Smart prep aorta at the level of the celiac artery with trigger set to 100 HU						
AXIAL ARTERIAL	Axial	MPR	Soft Tissue	2 mm	2 mm	Lung base through pubic symphysis
AXIAL DELAY 120 SECONDS	Axial	MPR	Soft Tissue	2 mm	2 mm	Lung base through pubic symphysis
COR MPR x 2 (Arterial and delay)	Coronal	MPR	Soft Tissue	2 mm	2 mm	
SAG MPR x 2 (Arterial and delay)	Sagittal	MPR	Soft Tissue	2 mm	2 mm	

Notes:

- Contrast: Per site, 100 cc.
- Craniocaudal scan direction.
- No oral contrast.