

MR Abdomen (Liver or Routine)

Last Updated: 12/2025

Name	Plane	Sequence
COR T2 Single Shot	Coronal	T2 Single Shot (AKA: SS TSE / HASTE)
AX T2 Single Shot	Axial	T2 Single Shot (AKA: SS TSE / HASTE)
AX T2 FATSAT	Axial	T2 FAT SAT FSE (AKA: TSE)
AX DWI	Axial	DWI (b=1000)
AX T1	Axial	DUAL ECHO IN/OUT PHASE GRE
AX DYNAMIC GAD	Axial	3D GRE FAT SAT LAVA MULTIPHASE (AKA: THRIVE / VIBE) (INCLUDE PRE-CONTRAST IN SAME RUN)
COR GAD	Coronal	LAVA DELAYED (AKA: THRIVE / VIBE)

Notes:

- The LAVA sequence is set on multiphase 6 sets. Run the first one before contrast injection.
- Pre, early arterial and/or late arterial, portal venous phase, equilibrium 3min, and 5 min delay.