

## MR Liver (Eovist)

Last Updated: 9/2/2025

Name	Plane	Sequence
<b>COR T2 Single Shot</b>	Coronal	T2 Single Shot (AKA: SS TSE / HASTE)
<b>AX T2 Single Shot</b>	Axial	T2 Single Shot (AKA: SS TSE / HASTE)
<b>AX T2 FATSAT</b>	Axial	T2 FAT SAT FSE (AKA: TSE)
<b>AX DWI</b>	Axial	DWI (b=1000)
<b>AX T1</b>	Axial	DUAL ECHO IN/OUT PHASE GRE
<b>Additional 10 and 20 minute delayed post contrast axial.</b>		
<b>AX DYNAMIC GAD</b>	Axial	3D GRE FAT SAT LAVA MULTIPHASE (AKA: THRIVE / VIBE) (INCLUDE PRE-CONTRAST IN SAME RUN)
<b>COR GAD</b>	Coronal	LAVA DELAYED (AKA: THRIVE / VIBE)

### Notes:

- Routine liver with additional 10 and 20 minute delayed post contrast axial.
- The LAVA sequence is set on multiphase 6 sets. Run the first one before contrast injection.
- Pre, early arterial and/or late arterial, portal venous phase, equilibrium 3min, and 5 min delay.