

MR Pancreas

Last Updated: 9/4/2025

Name	Plane	Sequence	Slice
AX T2 Single Shot	Axial	T2 Single Shot (AKA: SS TSE / HASTE)	
COR T2 Single Shot	Coronal	T2 Single Shot (AKA: SS TSE / HASTE)	
AX DWI	Axial	DWI (b=1000)	
AX T1	Axial	DUAL ECHO IN/OUT PHASE GRE	
AX T2 FATSAT	Axial	T2 FAT SAT FRFSE (AKA: RESTORE / DRIVE)	
COR MRCP	Coronal	3D FRFSE (AKA: RESTORE / DRIVE) RESP TRIGGER	
MRCP RADIAL	Radial	MRCP 9 SHOT HEAVY T2 WEIGHTED SSFSE (AKA: SS TSE / HASTE)	THIN AND/OR THICK SLAB
AX DYNAMIC GAD	Axial	LAVA 3D GRE FAT SAT Dynamic (INCLUDE PRE-CONTRAST) (AKA: THRIVE / VIBE) [up to 5 min]	
COR GAD	Coronal	LAVA DELAY (AKA: THRIVE / VIBE)	

Notes:

- The LAVA sequence is set on multiphase 6 sets. Run the first one before contrast injection.
- Pre, art, late art, portal venous phase, equilibrium 3min, and 5 min delay.