

MR Pelvis Perianal Fistula

Last Updated: 12/2025

Name	Plane	Thick	Space	Comments
SCOUT				
3D T2W TSE	Sagittal	3	-0.15	If site able to do a 3D T2
2D T2W TSE	Axial	4	0.4	
2D T2W SSTSE FS	Axial	4	0.4	
2D T1W IP/OP	Axial	5	1	
2D T2W TSE	Coronal	4	0.4	
3D T1W DIXON POST	Axial	3	-1.5	
3D T1W DIXON POST	Sagittal	3	-1.5	
OPTIONAL				
3D T1W DIXON POST	Coronal	3	-1.5	
2D DWI FS	Axial	5	0.5	b50, 500, 1000, ADC
CONTRAST				
TYPE				
DOSE				
DELAY				

Notes:

Recommend placement of gauze in gluteal cleft to separate near anus

Scan Range: Small FOV from mid rectum through anus for all series