

MR Urethra

Last Updated: 12/2025

Name	Plane	Sequence	Slice	Gap	FOV
AX T2 Single Shot	Axial	T2 Single Shot FSE (AKA: SS TSE / HASTE)	3-4 mm	10%	
COR T2 Single Shot	Coronal	T2 Single Shot FSE (AKA: SS TSE / HASTE)	3-4 mm	10%	
AX T2	Axial	T2 FSE (TSE)	3-4 mm	10%	Small FOV from bladder to perineum.
COR T2	Coronal	T2 FSE (TSE)	3-4 mm	10%	Small FOV centered on urethra.
SAG T2	Sagittal	T2 FSE (TSE)	3-4 mm	10%	Small FOV centered on urethra.
AX T1	Axial	DUAL ECHO IN/OUT PHASE	3-4 mm	10%	Small FOV
AX DYNAMIC GAD	Axial	LAVA DYNAMIC (AKA: THRIVE / VIBE) (INCLUDE PRE-CONTRAST)	3-4 mm	10%	Small FOV
SAG GAD	Sagittal	POST CONTRAST LAVA (AKA: THRIVE / VIBE)	3-4 mm	10%	Small FOV