

**Probes:** L12-3, eL18-4, C8-5 or appropriate transducer depending on patient size

**Image 2D, color and Doppler Waveform**

Arteries imaged: obtain PSV reading

- Innominate RT side only
- Subclavian
- Axillary
- Brachial
- Radial (if needed)
- Ulna (if needed)

On **normal** appearing Unilateral exams, image the Mid Subclavian Artery and Axillary A and Brachial A on the contralateral side to compare waveforms.

On **abnormal** appearing Unilateral exams image at and below the abnormal area on the contralateral side to compare waveforms of the abnormal side.

**Complete worksheet and scan into PACS**  
**Document Triphasic, Biphasic or Monophasic flow**  
**Document PSV**