

**Probes:** L12-3, eL18-4, C8-5 or appropriate transducer depending on patient size  
**image Longitudinal with Gray scale, Color, and Doppler Waveform**

Arteries imaged: Obtain PSV readings

- Common Iliac
- Internal Iliac if can be seen
- External Iliac
- Common Femoral
- Central Profunda Femoris artery
- Superficial Femoral – central, mid and peripheral
- Popliteal –central and peripheral
- Anterior Tibial
- Posterior Tibial –ACR requires only peripheral
- Peroneal
- Dorsalis Pedis (DPA)

On **normal** appearing Unilateral exams image the CFA, Mid FA and POP A on the contralateral side to compare waveforms and PSV's.

On **abnormal** appearing Unilateral exams image **at and below** the abnormal area on the contralateral side to compare waveforms and PSV's of the abnormal side.

If clot is seen CIA, please image into the AO for extension.

**Complete worksheet and scan into PACS**  
**Document Triphasic, Biphasic, or Monophasic flow**  
**Document PSV**